BADABULLE

FR Notice d'utilisation
NL Gebruiksaanwijzing
EN Instruction manual
IT Libretto d'utilizzazione
ES Instrucciones de uso
PT Instruções utilização
DE Anleitungen zur benutzung
CZ Uživatel
AR تعليمات التشغيل

Réf. B037200

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EN OPERATING INSTRUCTIONS

Before using, please disinfect the probe at first. To switch on, press the ON/OFF button next to the display; a short beep will sound. When the letters «Lo» and a flashing «°C» display, the thermometer is now ready for use. If the ambient temperature is below 32°C, then «Lo °C» will appear on the LCD and if it is more than 42.9°C, then «Hi °C» will appear on the LCD. During the reading, the current temperature is displayed continuously and the «°C» symbol flashes. The measurement is completed when a constant temperature value has been reached. So, a beep will sound ten times, and the «°C» symbol will stop flashing. The highest measured temperature appears on the LCD. When the measurement is completed, please switch the thermometer off by pressing the ON/OFF button. After the temperature has been displayed, the thermometer will shut off automatically after 10 minutes.

MEMORY FUNCTION

Turn on the thermometer; the last measured temperature will display on the LCD for approximately 2s after the self-check. The reading will be overwritten after the LCD shows «37.0°C», i.e. the thermometer is ready for a next measurement. Methods of measuring temperature It is important to remember that the body temperature reading depends on the site where it is measured. For this reason, the measurement site must always be specified in order to ensure that a correct temperature reading is recorded.

In the rectum (rectal)

This is the most accurate method from a medical point of view. Insert carefully the thermometer tip into the rectum for a maximum of 2 cm. The usual measuring time is approximately 9 to 11 seconds.

Under the arm (axillary)

The measurement can fluctuate by about 0.5°C to 1.5°C from rectal temperature readings. The usual measuring time for this method is approximately 23 to 29 seconds. It should be noted, however, that an exact reading cannot be obtained if, for example, the armpits have been allowed to cool. If this is the case, we recommend extending the measuring time by around 5 minutes in order to obtain the most precise information.

In the mouth (oral)

The oral temperature is 0.3°C to 0.8°C lower than the rectal temperature. To ensure that readings as accurate as possible, place the thermometer tip to the left or right of the root of the tongue. The thermometer tip must have constant contact with the tissue, keep the mouth closed. The usual measuring time is approximately 15 to 19 seconds. We advise you to extend the measuring time by 3 minutes after the beep.

CLEANING AND DISINFECTION

Apply a disinfectant with a damp cloth

SAFETY PRECAUTIONS

Do not allow the device to come into contact with hot water.
Do not expose to high temperatures or direct sunlight.
Do not drop the thermometer. It is neither shock-proof nor impact-resistant.
Do not modify this device without the authorization of the manufacturer.
Do not bend or open the device (except the battery compartment).
Do not clean with thinners, petrol or benzene. Only clean with water or disinfectant.
The thermometer contains small parts (battery, battery compartment) which can be swallowed by children. For this reason, do not leave the thermometer unattended in the hands of children.
Avoid bending the tip by more than 45 degrees.
If the ambient temperature is over 35 °C, dip the thermometer tip in cold water for approx. 5 to 10 seconds prior to measuring the temperature.
Persistent fever, in particular in children, has to be treated by a doctor - please get in touch with your doctor!
Do not use near strong electromagnetic fields, i.e. Keep it away from any radio systems and mobile phones.

BATTERY REPLACEMENT

The battery needs replacing when the «» or «» battery symbol appears on the right of the LCD. Remove the battery cover and replace it with a battery (preferably non-mercury) of the same type. Please note: The «+» sign up and «-» sign down.

PRODUCT DISPOSAL

Please ensure environmental protection. Batteries do Not belong in the domestic waste. Please hand them in at collection point. The alkaline battery or fuel cell may lead to excessive temperatures, fire or explosion.

This symbol on products and/or accompanying documents means that consumed electronic products must not be mixed with conventional domestic waste. Take these products to the corresponding collection points for correct treatment and recycling, where they will be accepted free of charge. For more information on the closest collection point, please enquire with your local authorities.

TECHNICAL DATA

Type: maximum thermometer
Measurement range: 32.0-42.9°C
Measurement accuracy: +/-0.1°C
(35.5°C-95.9°C)
Storage/transportation temperature: -25-55°C, <=95%RH
Ambient temperature during use: 5-35°C, <=80%RH
Battery type: Alkaline battery, type LR41, 1.5V, service life minimum 100 hours under continuous operation.
Weight: approx. 13g
Atmospheric pressure: 700-1060hPa
Shell life: 3 years

LEGAL REQUIREMENTS AND GUIDELINES

This product complies with the European Directive for Medical Devices 93/42/EEC and carries the CE mark. The CE marking confirms that this is a medical device with a measuring function in the sense of the Medical Devices Act which has undergone a conformity assessment procedure. A Notified Body confirms that this product fulfills all the appropriate statutory regulations.

WARRANTY

This product is warranty for 1 year from the date of leave factory. Damage resulting from incorrect use or abuse is not covered by the warranty. Battery and packaging are excluded from the warranty. Claims beyond this, including claims for damages, are excluded. If you find that the thermometer is defective and not in good function, please firstly check the battery before sending in for repair.