# BABY MOOV

Réfs. : A032002 - A032004 - A032007 - A032008

- NOTICE D'UTILISATION DU THERMOMÈTRE RECTAL NUMERIQUE ET DU MOUCHE BÉBÉ
- INSTRUCTIONS MANUAL FOR DIGITAL RECTAL THERMOMETER AND NASAL ASPIRATOR
- GEBRAUCHSANWEISUNG FÜR DAS DIGITALE REKTALTHERMOMETER UND DEN NASENSAUGER
- GEbruIksaanwIzIng DIGITale RECTale THERMOMETER EN DE BABYSnUITER
- INSTRUCCIONES DE USO TERMÓMETRO RECTAL DIGITAL Y DEL SACAMOCOS
- ISTRUZIONI PER L'USO DEL TERMOMETRO RETTALE DIGITALE E DELL'ASPIRATORE NASALE
- INSTRUÇÕES DE UTILIZAÇÃO DO TERMÓMETRO RECTAL DIGITAL E DO ASPIRADOR NASAL PARA BEBÉ
- NÁVOD K POUŽITÍ DIGITÁLNÍHO REKTÁLNÍHO TEPLOMĚRU A NOSNÍ ODSÁVAČKY
- A DIGITÁLIS REKTÁLIS HŐMÉRŐ ÉS AZ ORRSZÍVÓ HASZNÁLATÁRA VONATKOZÓ UTASÍTÁSOK
- INSTRUCȚIUNI DE UTILIZARE TERMOMETRU RECTAL DIGITAL ȘI A ASPIRATORULUI NAZAL PENTRU BEBELUȘI
- INSTRUKCJA OBSŁUGI CYFROWEGO TERMOMETRU DOODBYTNICZEGO I ASPIRATORA DLA NIEMOWIĄT

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EN.

## IMPORTANT. READ INSTRUCTIONS BEFORE USE AND KEEP FOR FUTURE REFERENCE

FEATURES	
- Temperature Range: 32 °C - 42,9 °C (90 °F - 109,9 °F)	- Battery life: about 250 hours or 1000 times.
- Accuracy: ± 0,1 °C (35 °C - 42 °C) ; ± 0,2 °C (< 35 °C)	- Flexible tip for comfort and safety.
- Digital LCD display.	- Beeps when reaches highest temperature.
- Battery: 1,5 VDC (LR41) alkaline.	- Auto-off in approximately 10 minutes.
- Power consumption: 0,15 mW.	- Last measurement memory.
- Dimensions: 118 mm (L)x 25 mm (W) x 10 mm (H).	- Low battery indicator ▼▼.
	- Weight with battery: 10 g.

<b>IP22</b> Dust and water resistance	 Operating instructions	 Type BF
---------------------------------------	--	---

#### HOW TO TAKE TEMPERATURE?

- Press power switch, the screen shows the full segment 188.8 °E with a sound of "B!", which indicates sufficient power and normal self-test routine as well as for measurement.
- Memory function: after the full segment, the last measured value with the letter "M" are displayed, then to show 37 °C / 98,6 °F.
- Release the button, then displays « Lo °C » or « Lo °F » and ready to measure. During measurement the symbol « °C » or « °F » will flash while the temperature is diming. This will also happen when temperature increase remains below 0,1 °C / 0,2 °F for about 16 seconds. If temperature remains below 32 °C / 90 °F, the symbol « °C » or « °F » will flash for about 60 seconds. Then it will auto-off for saving power.
- A buzzing tone will be heard for about 10 seconds which means the finish of measurement. Normal alarm: sound of « .. bip .. bip .. bip .. bip ». « .. bip .. bip .. bip » means normal temperature that below 37,8 °C / 100 °F. Fever alarm: sound of « ..bip..bip..bip..bip » means fever that the temperature is above 37,8 °C / 100 °F.
- During measurement: LCD will display « Lo °C » or « Lo °F » if the measured temperature is below 32,0 °C (90 °F), or will display « H °C » or « H °F » if the measures temperature is above 43 °C (109,9 °F).
- Turn off the product after slightly pressing power switch, or leave it for 9 minutes, it will auto-off for saving power.
- Children should be under close supervision at all times while using this product. After monitoring the temperature, the unit should be removed from the child immediately.
- Normal temperature: about 37,8 °C / 99,7 °F.
- °C » or « °F » exchange function: if the power switch was pushed more than 2 sec after power-on, the unit °C (°F) will be exchanged.

#### DIFFERENT METHODS OF USE

- You are advised to measure your temperature at specific time in morning, if it is earlier or later, the temperature would be lower or higher.
- Used in mouth.** Please close your mouth for 1 minute for stability of temperature, breath calmly through nose. Respiratory air would affect the accuracy of measure of your mouth. Normal temperature: about 37 °C / 98,6 °F.
- Used in anus of rectum.** This is the way that can get the most steady temperature and suits for both infants and children. The sensor shall be inserted for 2-3 cm. Vaseline can be used for easier inserting.
- Normal temperature: about 37,8 °C / 99,7 °F.
- Used under armpit.** Clean the thermometer with a dry cloth and put it under armpit. Press and fix it with the arm. The temperature value of this measurement is not so accurate that not recommended for accurate measurement.
- The scrap of product and battery are to be disposed according to relevant national regulations.
- Notes: the normal temperatures are different from the different human bodies and measuring time. Sterilization is required before every measurement.

#### BATTERY REPLACEMENT

Battery must be replaced if ▼▼ symbol appears during measurement. Make sure a battery of the same type is on hand and proceed as follows: Remove battery compartment lid and change with a button battery (LR41). Ensure battery is correctly positioned with "+" facing up.
**WARNING** keep battery out of children's reach. Ingestion can cause chemical burns, perforation of soft tissues and death. Serious burns can occur within two hours of ingestion. Seek medical attention immediately. Battery should not be charged or placed into extreme heat as it may explode.

#### CLEANING AND MAINTENANCE

- The thermometer is waterproof and can be cleaned by water or sterilization lotion with soft cloth. Sterilization lotion below is recommended:
**Note**: Density - Medical alcohol: 95 %
**Notes**: other lotion, density or cleaning method may damage the thermometer.
**Baby nasal aspirator**: Clean regularly using warm soapy water only, ensuring it is thoroughly rinsed. Dry immediately. Do not use the product if any parts are broken or missing.
Any serious accident arising from the use of this device should be notified to the manufacturer and competent authority in the member state where the user lives.

#### CAUTIONS

- How to use the baby nasal aspirator**: can be used from birth. Unblocks the nose. Keep out of reach of children. Do not expose to high temperatures or direct sunlight.
- Store the unit in the protective case while not in use.
- Do not store or keep in direct sunlight, or in high temperatures.
- Keep away from strong electromagnetic source to avoid interference.
- Do not open the thermometer except to replace battery.
- Do not drop the unit, as this may ruin the unit.
- It shall be operated in an environment at the temperature of 10 °C / 50 °F and the humidity of 40-80 % while stored and transported in an environment at the temperature of -10 °C - 40 °C (14 °F - 104 °F) and humidity of 30 - 90 %.
- Never put it in boiled water for sterilizing.
- Please take out the battery if not use the unit for long time.

**REMARKS** :
It is suitable for measuring temperature of human body and necessary of families.
It is designed and manufactured complying with relevant requirements in ISO 80601-2-56.
It has passed EMC emission test, immunity test according to EN60601-1-2 and biocompatibility test according to ISO 10993-1.
According to the degree of protection against electric shock, it is classified to BF, in accordance with EN60601-1.
Atmospheric pressure range is 70–106Pa during operation
Atmospheric pressure range is 70–106Pa during transportation
**WARNING**: Modify/repair this device is not allowed.
**WARNING**: do not modify this equipment without authorization of the manufacturer.
**WARNING**: this equipment is modified, appropriate inspection and testing must be.
Please replace battery before claim for replacing a new product.
**Caution** Be sure that both the battery and the thermometer are kept away from children.
Unwanted electrical and electronic items are collected separately. Do not discard unwanted electrical and electronic items with unsorted municipal waste, but ensure that they are collected separately.
These instructions were revised in april 2024.

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